


DEDICATED CONTRAST MEDIA INJECTOR FOR MAMMOGRAPHY **APO 100M**

*The 1st Contrast Media Injector,
dedicated to Mammography rooms!*



APOLLO RT Technology
Service People

Compact and Easy to Move :

APO100M can be easily moved close to the patient when needed, and placed aside during standard exams.

Slow Speed Keys

Slowly moves forward & backward injector plunger. Can be used for air purging, or confirming needle placement.

Built in heat maintainer

Built in heat maintainer keeps **contrast media** at convenient temperature



Remote Hand Switch:

Allows the medical staff to be close to the patient, and monitor vein patency during whole injection

Waterproof Keys

All keys are waterproof designed with soft touch.

Safety Confirmation Key

APO100M injector head has a built-in Safety Confirmation Key. Injection can only be started after proper checking is performed, and after this key is pressed. This is preventing:

- ♦ Air injection;
- ♦ Undesired injections, in case Start Switch is activated;


Integrated AC Power Supply to replace Battery:

good for the environment, efficient, less maintenance, less cost, etc. compared to using battery.

- ♦ Better for environment;
- ♦ More efficient;
- ♦ Less costly;
- ♦ No maintenance


Reliable Casters:

Large stand wheels with locking capabilities



User *Friendly* Console

- Large & Bright LCD color touch screen, with wide vision angle. Parameters and injection status are clearly displayed.
- User friendly interface, so that any access can be handled by one hand.
- Large font size numeric keypad is convenient for data entering.



Real Time Monitoring of pressure graph, flow rate & injected volume

- Real time pressure is displayed during injection.
- If needed pressure to reach Set Flow Rate is exceeding set Pressure Limit, APO100M automatically decreases Flow Rate.
- If Actual Pressure is higher than set Pressure limit, APO100M automatically stops.
- If Flow Rate or Injected Volume are higher than set values, APO100M stops immediately.



Multi Injection modes are fulfilling various needs

- Two Modes are available: single phase injection, dual phase injection.
- UNIQUE Technology: **Body Weight Protocol Injection.** By entering the patient body weight and contrast iodine concentration, APO100M will automatically calculate the optimal injection protocol. This is ensuring simple manipulation process, and maximizing work efficiency, while improving injection safety by preventing wrong parameter settings.


APO100M disposables are made up of highly clear material, with high pressure resistance.

All production steps of APO100M disposables are controlled by Scientific Quality Management.

This is, indeed, the Injector YOU NEED in your mammographic room.

Technical Data

Flow Rate	0.1-5.0ml/s (0.1ml/s increment)
Volume	1-200ml (1 ml increment)
Upper Pressure Limit	300PSI
Injection Mode	Multi phase
Heating	Built-in Temperature Control Heating(38° +/- 3°)
Rated Voltage	220AV,50Hz/110AV,60Hz
Power Output Or Consumption	100 VA
Installation	Pedestal(standard)
Consumables	100ml syringe/200ml syringe



200ml Syringe



100ml Adaptor



100ml Syringe

Caution: Safety & Performance of APO100M cannot be warranted if used with syringes or disposables other than approved by Company. Specifications & Appearance, as shown in this brochure, can be changed by Company without notice.

APOLLO RT Co., Ltd.

Unit 1101, 11/F, Tower 1, Cheung Sha Wan Plaza, 833 Cheung Sha Wan Road, Lai Chi Kok, Kowloon, Hong Kong
Tel: 852-2781-3481
e-mail: contact@apollo-rt.com

Manufactured by

APOSTAR Medical (Shanghai) Limited

Add: Suit B, 2/F, Building 3, No.1151 Lianxi Road, Pudong 201204 Shanghai